



63 East 1280 North Tooele UT 84074  
dependableoxygen.com | (435)882-0055

## Patient Cleaning Procedures-Humidifier Bottle

### Daily

- Clean humidifier bottle using warm soapy water. (You can also place on top shelf of dish washer.)
- Rinse with hot water.
- Refill with distilled water only, don't exceed the maximum line.

### Disinfect the bottle 1-2 times per week.

- Soak in a solution of 1 part white vinegar, 3 part water, and let soak for 30 minutes.
- Rinse with hot water.
- Refill with Distilled water.

I have been instructed how to install the oxygen humidifier, and understand how the cleaning procedures.

---

Patients Signature

---

Date

---

D.O.C Signature

---

Date